

Varenavn:

**ASAHI wires**









Produktgruppe:

Nevrovaskulære Guide Wires

Produsent:

ASAHI



Produkt	Varenummer (code)	Produktbeskrivelse	Wire lengde (cm)	Hydrofil Coating	Polymer sleeve	Spring Coil lengde (cm)	Tipp form	Røntgentett lengde (cm)
CHIKAI 008	WAIN-CKI-008-200	Asahi Chikai 008, 200 cm x 0.008"/0.010"	200	180	-	9	rett	9
CHIKAI 10	WAIN-CKI-10-200	Asahi Chikai 10, 200 cm x 0.010"	200	170	-	9,5	rett	3
	WAIN-CKI-10-300	Asahi Chikai 10, 300 cm x 0.010"	300	170	-	9,5	rett	3
CHIKAI X 010	AIN-CKX-10-200-R	Asahi Chikai X 010, 200 x 0.010" med Re-shapeable 25° bøy	200	180	150	16		16
CHIKAI	WAIN -CKI-200	Asahi Chikai, 200 cm x 0.014" rett	200	170	-	30	rett	5
	WAIN -CKI-200-RC	Asahi Chikai, 200 cm x 0.014" med bøy	200	170	-	30		5
	WAIN-CKI-300	Asahi Chikai, 300 cm x 0.014" rett	300	170	-	30	rett	5
	WAIN-CKI-300-RC	Asahi Chikai, 300 cm x 0.014" med bøy	300	170	-	30		5
CHIKAI Black	WAIN-CKI-200-BS	Asahi Chikai Black 200 cm x 0,014" med bøy	200	170	150	30		3
	WAIN-CKI-200 BA	Asahi Chikai Black 200 cm x 0,014" med 90° bøy	200	170	150	30		3
CHIKAI Black 14 soft tip	AIN-CKI-200-B-SFT	Asahi SION black, J-bøy, 180 cm x 0.014" med bøy, myk tupp	200	170	150	30		3
CHIKAI Black 18	WAIN-CKI-18-200-BS	Asahi SION black, J-bøy, 300 cm x 0.018" med bøy	200	170	150	34		5
CHIKAI Black 18 soft tip	AIN-CKI-18-200-SFT	Asahi SION black, J-bøy, 300 cm x 0.018" med bøy, myk tupp	200	170	150	34		5



Act One Technology: enhances tip durability and reduces whipping of the wire due to superior torque control.